



SLC14A1 Antibody, Biotin conjugated

Product Code	CSB-PA615678LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13336
Immunogen	Recombinant Human Urea transporter 1 protein (1-51AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Urea transporter 1 (Solute carrier family 14 member 1) (Urea transporter, erythrocyte), SLC14A1, HUT11 JK RACH1 UT1 UTE
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC14A1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.