



LAMC2 Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA618803LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q13753 |
| Immunogen | Recombinant Human Laminin subunit gamma-2 protein (417-588AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Laminin subunit gamma-2 (Cell-scattering factor 140 kDa subunit) (CSF 140 kDa subunit) (Epiligrin subunit gamma) (Kalinin subunit gamma) (Kalinin/nicein/epiligrin 100 kDa subunit) (Ladsin 140 kDa subunit) (Laminin B2t chain) (Laminin-5 subunit gamma) (Large adhesive scatter factor 140 kDa subunit) (Nicein subunit gamma), LAMC2, LAMB2T LAMNB2 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | LAMC2 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |