



LAMC2 Antibody, Biotin conjugated

Product Code	CSB-PA618803LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q13753
Immunogen	Recombinant Human Laminin subunit gamma-2 protein (417-588AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Laminin subunit gamma-2 (Cell-scattering factor 140 kDa subunit) (CSF 140 kDa subunit) (Epiligrin subunit gamma) (Kalinin subunit gamma) (Kalinin/nicein/epiligrin 100 kDa subunit) (Ladsin 140 kDa subunit) (Laminin B2t chain) (Laminin-5 subunit gamma) (Large adhesive scatter factor 140 kDa subunit) (Nicein subunit gamma), LAMC2, LAMB2T LAMNB2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LAMC2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.