



TNIP1 Antibody, HRP conjugated

Product Code	CSB-PA619959LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15025
Immunogen	Recombinant Human TNFAIP3-interacting protein 1 protein (526-636AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	TNFAIP3-interacting protein 1 (A20-binding inhibitor of NF-kappa-B activation 1) (ABIN-1) (HIV-1 Nef-interacting protein) (Nef-associated factor 1) (Naf1) (Nip40-1) (Virion-associated nuclear shuttling protein) (VAN) (hVAN), TNIP1, KIAA0113 NAF1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	TNIP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.