



PAFAH1B3 Antibody, FITC conjugated

Product Code	CSB-PA619963LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15102
Immunogen	Recombinant Human Platelet-activating factor acetylhydrolase IB subunit gamma protein (103-187AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Platelet-activating factor acetylhydrolase IB subunit gamma (EC 3.1.1.47) (PAF acetylhydrolase 29 kDa subunit) (PAF-AH 29 kDa subunit) (PAF-AH subunit gamma) (PAFAH subunit gamma), PAFAH1B3, PAFAHG
Immunogen Species	Homo sapiens (Human)
Research Area	Cardiovascular
Target Names	PAFAH1B3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.