



SIM2 Antibody, HRP conjugated

Product Code	CSB-PA621763LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q14190
Immunogen	Recombinant Human Single-minded homolog 2 protein (347-433AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Single-minded homolog 2 (Class E basic helix-loop-helix protein 15) (bHLHe15), SIM2, BHLHE15
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	SIM2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.