



USF2 Antibody, Biotin conjugated

Product Code	CSB-PA621888LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q15853
Immunogen	Recombinant Human Upstream stimulatory factor 2 protein (21-162AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Upstream stimulatory factor 2 (Class B basic helix-loop-helix protein 12) (bHLHb12) (FOS-interacting protein) (FIP) (Major late transcription factor 2) (Upstream transcription factor 2), USF2, BHLHB12
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	USF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.