



NSUN2 Antibody, HRP conjugated

Product Code	CSB-PA626626LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08J23
Immunogen	Recombinant Human tRNA (cytosine(34)-C(5))-methyltransferase protein (432-509AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	tRNA (cytosine(34)-C(5))-methyltransferase (EC 2.1.1.203) (Myc-induced SUN domain-containing protein) (Misu) (NOL1/NOP2/Sun domain family member 2) (Substrate of AIM1/Aurora kinase B) (tRNA (cytosine-5-)-methyltransferase) (tRNA methyltransferase 4 homolog) (hTrm4), NSUN2, SAKI TRM4
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NSUN2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.