



TRIM71 Antibody, HRP conjugated

Product Code	CSB-PA652403LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q2Q1W2
Immunogen	Recombinant Human E3 ubiquitin-protein ligase TRIM71 protein (125-209AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase TRIM71 (EC 2.3.2.27) (Protein lin-41 homolog) (RING-type E3 ubiquitin transferase TRIM71) (Tripartite motif-containing protein 71), TRIM71, LIN41
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TRIM71
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.