




ZAR1 Antibody

Product Code	CSB-PA653624LA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ZU46
Immunogen	Recombinant Arabidopsis thaliana Receptor protein kinase-like protein ZAR1 protein (167-278AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Receptor protein kinase-like protein ZAR1 (EC 2.7.11.1) (Protein ZYGOTIC ARREST 1), ZAR1
Immunogen Species	Arabidopsis thaliana
Research Area	Others
Target Names	ZAR1

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: ZAR1 antibody at 5.32µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 42 kDa

Observed band size: 42 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.