



ZAR1 Antibody, FITC conjugated

Product Code	CSB-PA653624LC01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9ZU46
Immunogen	Recombinant Arabidopsis thaliana Receptor protein kinase-like protein ZAR1 protein (167-278AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Receptor protein kinase-like protein ZAR1 (EC 2.7.11.1) (Protein ZYGOTIC ARREST 1), ZAR1
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Research Area	Others
Target Names	ZAR1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.