



# WLS Antibody, HRP conjugated

|                            |                                                                                                                                                                            |
|----------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA716359LB01HU                                                                                                                                                         |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                              |
| <b>Uniprot No.</b>         | Q5T9L3                                                                                                                                                                     |
| <b>Immunogen</b>           | Recombinant Human Protein wntless homolog protein (123-232AA)                                                                                                              |
| <b>Raised In</b>           | Rabbit                                                                                                                                                                     |
| <b>Species Reactivity</b>  | Human                                                                                                                                                                      |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                      |
| <b>Form</b>                | Liquid                                                                                                                                                                     |
| <b>Conjugate</b>           | HRP                                                                                                                                                                        |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                                                                                           |
| <b>Purification Method</b> | >95%, Protein G purified                                                                                                                                                   |
| <b>Isotype</b>             | IgG                                                                                                                                                                        |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                 |
| <b>Alias</b>               | Protein wntless homolog (Integral membrane protein GPR177) (Protein evenness interrupted homolog) (EVI) (Putative NF-kappa-B-activating protein 373), WLS, C1orf139 GPR177 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                                                                       |
| <b>Research Area</b>       | Signal Transduction                                                                                                                                                        |
| <b>Target Names</b>        | WLS                                                                                                                                                                        |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                                                                                                |