



WLS Antibody, Biotin conjugated

Product Code	CSB-PA716359LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5T9L3
Immunogen	Recombinant Human Protein wntless homolog protein (123-232AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein wntless homolog (Integral membrane protein GPR177) (Protein evenness interrupted homolog) (EVI) (Putative NF-kappa-B-activating protein 373), WLS, C1orf139 GPR177
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	WLS
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.