



# RNF213 Antibody, HRP conjugated

<b>Product Code</b>	CSB-PA717453LB01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q63HN8
<b>Immunogen</b>	Recombinant Human E3 ubiquitin-protein ligase RNF213 protein (202-350AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	HRP
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	E3 ubiquitin-protein ligase RNF213 (EC 2.3.2.27) (EC 3.6.4.-) (ALK lymphoma oligomerization partner on chromosome 17) (Mysterin) (RING finger protein 213) (RING-type E3 ubiquitin transferase RNF213), RNF213, ALO17 C17orf27 KIAA1554 KIAA1618 MYSTR
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Others
<b>Target Names</b>	RNF213
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.