



RNF213 Antibody, Biotin conjugated

Product Code	CSB-PA717453LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q63HN8
Immunogen	Recombinant Human E3 ubiquitin-protein ligase RNF213 protein (202-350AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase RNF213 (EC 2.3.2.27) (EC 3.6.4.-) (ALK lymphoma oligomerization partner on chromosome 17) (Mysterin) (RING finger protein 213) (RING-type E3 ubiquitin transferase RNF213), RNF213, ALO17 C17orf27 KIAA1554 KIAA1618 MYSTR
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RNF213
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.