



CEP135 Antibody, FITC conjugated

Product Code	CSB-PA721063LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q66GS9
Immunogen	Recombinant Human Centrosomal protein of 135 kDa protein (94-195AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Centrosomal protein of 135 kDa (Cep135) (Centrosomal protein 4), CEP135, CEP4 KIAA0635
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	CEP135
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.