



FREM2 Antibody, Biotin conjugated

Product Code	CSB-PA722558LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5SZK8
Immunogen	Recombinant Human FRAS1-related extracellular matrix protein 2 protein (1986-2127AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	FRAS1-related extracellular matrix protein 2 (ECM3 homolog), FREM2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	FREM2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.