



ATAT1 Antibody, HRP conjugated

Product Code	CSB-PA725542LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q5SQI0
Immunogen	Recombinant Human Alpha-tubulin N-acetyltransferase 1 protein (194-238AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Alpha-tubulin N-acetyltransferase 1 (Alpha-TAT) (Alpha-TAT1) (TAT) (EC 2.3.1.108) (Acetyltransferase mec-17 homolog), ATAT1, C6orf134 MEC17
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	ATAT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.