



SPTBN1 Antibody, FITC conjugated

Product Code	CSB-PA733819LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q01082
Immunogen	Recombinant Human Spectrin beta chain, non-erythrocytic 1 protein (2096-2256AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Spectrin beta chain, non-erythrocytic 1 (Beta-II spectrin) (Fodrin beta chain) (Spectrin, non-erythroid beta chain 1), SPTBN1, SPTB2
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	SPTBN1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.