



RHOU Antibody, Biotin conjugated

Product Code	CSB-PA745331LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q7L0Q8
Immunogen	Recombinant Human Rho-related GTP-binding protein RhoU protein (1-48AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rho-related GTP-binding protein RhoU (CDC42-like GTPase 1) (GTP-binding protein-like 1) (Rho GTPase-like protein ARHU) (Ryu GTPase) (Wnt-1 responsive Cdc42 homolog 1) (WRCH-1), RHO, ARHU CDC42L1 G28K WRCH1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	RHO
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.