



PIK3AP1 Antibody, HRP conjugated

Product Code	CSB-PA747806LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6ZUJ8
Immunogen	Recombinant Human Phosphoinositide 3-kinase adapter protein 1 protein (293-404AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Phosphoinositide 3-kinase adapter protein 1 (B-cell adapter for phosphoinositide 3-kinase) (B-cell phosphoinositide 3-kinase adapter protein 1), PIK3AP1, BCAP
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	PIK3AP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.