



SCN10A Antibody, Biotin conjugated

Product Code	CSB-PA750910LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y5Y9
Immunogen	Recombinant Human Sodium channel protein type 10 subunit alpha protein (992-1099AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium channel protein type 10 subunit alpha (Peripheral nerve sodium channel 3) (PN3) (hPN3) (Sodium channel protein type X subunit alpha) (Voltage-gated sodium channel subunit alpha Nav1.8), SCN10A
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SCN10A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.