



CD300LG Antibody, Biotin conjugated

Product Code	CSB-PA751016LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6UXG3
Immunogen	Recombinant Human CMRF35-like molecule 9 protein (128-241AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CMRF35-like molecule 9 (CLM-9) (CD300 antigen-like family member G) (Triggering receptor expressed on myeloid cells 4) (TREM-4) (CD antigen CD300g), CD300LG, CLM9 TREM4
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD300LG