



PDXDC1 Antibody, FITC conjugated

Product Code	CSB-PA761311LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6P996
Immunogen	Recombinant Human Pyridoxal-dependent decarboxylase domain-containing protein 1 protein (681-780AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Pyridoxal-dependent decarboxylase domain-containing protein 1 (EC 4.1.1.-), PDXDC1, KIAA0251
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PDXDC1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.