



RBMS3 Antibody, FITC conjugated

Product Code	CSB-PA761548LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q6XE24
Immunogen	Recombinant Human RNA-binding motif, single-stranded-interacting protein 3 protein (372-437AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding motif, single-stranded-interacting protein 3, RBMS3
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	RBMS3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.