



NLRP10 Antibody, HRP conjugated

Product Code	CSB-PA773044LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86W26
Immunogen	Recombinant Human NACHT, LRR and PYD domains-containing protein 10 protein (589-655AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NACHT, LRR and PYD domains-containing protein 10 (Nucleotide-binding oligomerization domain protein 8), NLRP10, NALP10 NOD8 PYNOD
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	NLRP10