



PDSS2 Antibody, FITC conjugated

Product Code	CSB-PA773060LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86YH6
Immunogen	Recombinant Human Decaprenyl-diphosphate synthase subunit 2 protein (292-387AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Decaprenyl-diphosphate synthase subunit 2 (EC 2.5.1.91) (All-trans-decaprenyl-diphosphate synthase subunit 2) (Candidate tumor suppressor protein) (Decaprenyl pyrophosphate synthase subunit 2), PDSS2, C6orf210 DLP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PDSS2