



# PDSS2 Antibody, Biotin conjugated

<b>Product Code</b>	CSB-PA773060LD01HU
<b>Storage</b>	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
<b>Uniprot No.</b>	Q86YH6
<b>Immunogen</b>	Recombinant Human Decaprenyl-diphosphate synthase subunit 2 protein (292-387AA)
<b>Raised In</b>	Rabbit
<b>Species Reactivity</b>	Human
<b>Tested Applications</b>	ELISA
<b>Form</b>	Liquid
<b>Conjugate</b>	Biotin
<b>Storage Buffer</b>	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
<b>Purification Method</b>	>95%, Protein G purified
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Alias</b>	Decaprenyl-diphosphate synthase subunit 2 (EC 2.5.1.91) (All-trans-decaprenyl-diphosphate synthase subunit 2) (Candidate tumor suppressor protein) (Decaprenyl pyrophosphate synthase subunit 2), PDSS2, C6orf210 DLP1
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Research Area</b>	Signal Transduction
<b>Target Names</b>	PDSS2