



PDSS2 Antibody, Biotin conjugated

Product Code	CSB-PA773060LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86YH6
Immunogen	Recombinant Human Decaprenyl-diphosphate synthase subunit 2 protein (292-387AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Decaprenyl-diphosphate synthase subunit 2 (EC 2.5.1.91) (All-trans-decaprenyl-diphosphate synthase subunit 2) (Candidate tumor suppressor protein) (Decaprenyl pyrophosphate synthase subunit 2), PDSS2, C6orf210 DLP1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PDSS2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.