



ABCA12 Antibody, HRP conjugated

Product Code	CSB-PA774804LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86UK0
Immunogen	Recombinant Human ATP-binding cassette sub-family A member 12 protein (45-158AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ATP-binding cassette sub-family A member 12 (ATP-binding cassette transporter 12) (ATP-binding cassette 12), ABCA12, ABC12
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	ABCA12
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.