




SERPINA11 Antibody

Product Code	CSB-PA801228LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86U17
Immunogen	Recombinant Human Serpin A11 protein (279-394AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Serpin A11, SERPINA11
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SERPINA11

Image



Western Blot

Positive WB detected in: MCF-7 whole cell lysate, Raji whole cell lysate, HepG2 whole cell lysate

All lanes: SERPINA11 antibody at 5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 47 kDa

Observed band size: 47 kDa