



LSR Antibody, HRP conjugated

Product Code	CSB-PA801252LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86X29
Immunogen	Recombinant Human Lipolysis-stimulated lipoprotein receptor protein (112-233AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Lipolysis-stimulated lipoprotein receptor, LSR, LISCH
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	LSR
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.