



GATAD2A Antibody, Biotin conjugated

Product Code	CSB-PA803156LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86YP4
Immunogen	Recombinant Human Transcriptional repressor p66-alpha protein (170-241AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transcriptional repressor p66-alpha (Hp66alpha) (GATA zinc finger domain-containing protein 2A), GATAD2A
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	GATAD2A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.