



HIPK1 Antibody, HRP conjugated

Product Code	CSB-PA803159LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86Z02
Immunogen	Recombinant Human Homeodomain-interacting protein kinase 1 protein (1007-1151AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeodomain-interacting protein kinase 1 (EC 2.7.11.1) (Nuclear body-associated kinase 2), HIPK1, KIAA0630 MYAK NBAK2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HIPK1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.