



# HIPK1 Antibody, Biotin conjugated

Product Code	CSB-PA803159LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q86Z02
Immunogen	Recombinant Human Homeodomain-interacting protein kinase 1 protein (1007-1151AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeodomain-interacting protein kinase 1 (EC 2.7.11.1) (Nuclear body-associated kinase 2), HIPK1, KIAA0630 MYAK NBAK2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	HIPK1