



CACNB2 Antibody, Biotin conjugated

Product Code	CSB-PA804482LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q08289
Immunogen	Recombinant Human Voltage-dependent L-type calcium channel subunit beta-2 protein (490-637AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Voltage-dependent L-type calcium channel subunit beta-2 (CAB2) (Calcium channel voltage-dependent subunit beta 2) (Lambert-Eaton myasthenic syndrome antigen B) (MYSB), CACNB2, CACNLB2 MYSB
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	CACNB2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.