



Wdr37 Antibody, HRP conjugated

Product Code	CSB-PA805557LB01MO
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8CBE3
Immunogen	Recombinant Mouse WD repeat-containing protein 37 protein (1-496AA)
Raised In	Rabbit
Species Reactivity	Mouse
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	WD repeat-containing protein 37, Wdr37, Kiaa0982
Immunogen Species	Mus musculus (Mouse)
Research Area	Others
Target Names	Wdr37
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.