



PMRT15 Antibody, Biotin conjugated

Product Code	CSB-PA812747LD01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8GWT4
Immunogen	Recombinant Arabidopsis thaliana Protein arginine N-methyltransferase 1.5 protein (513-637AA)
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein arginine N-methyltransferase 1.5 (AtPMRT15) (AtPMRT5) (EC 2.1.1.319) (Shk1 kinase-binding protein 1 homolog), PMRT15, PMRT5 SKB1
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Research Area	Others
Target Names	PMRT15
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.