



RBM45 Antibody, HRP conjugated

Product Code	CSB-PA814203LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IUH3
Immunogen	Recombinant Human RNA-binding protein 45 protein (170-244AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding protein 45 (Developmentally-regulated RNA-binding protein 1) (RB-1) (RNA-binding motif protein 45), RBM45, DRB1 DRBP1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RBM45
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.