



RBM45 Antibody, HRP conjugated

Product Code	CSB-PA814203LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IUH3
Immunogen	Recombinant Human RNA-binding protein 45 protein (170-244AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	RNA-binding protein 45 (Developmentally-regulated RNA-binding protein 1) (RB-1) (RNA-binding motif protein 45), RBM45, DRB1 DRBP1
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	RBM45