



SKA3 Antibody, Biotin conjugated

Product Code	CSB-PA816899LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8IX90
Immunogen	Recombinant Human Spindle and kinetochore-associated protein 3 protein (101-195AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Spindle and kinetochore-associated protein 3, SKA3, C13orf3 RAMA1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	SKA3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.