



DTNBP1 Antibody, Biotin conjugated

Product Code	CSB-PA822207LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96EV8
Immunogen	Recombinant Human Dysbindin protein (258-351AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Dysbindin (Biogenesis of lysosome-related organelles complex 1 subunit 8) (BLOC-1 subunit 8) (Dysbindin-1) (Dystrobrevin-binding protein 1) (Hermansky-Pudlak syndrome 7 protein) (HPS7 protein), DTNBP1, BLOC1S8
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	DTNBP1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.