



MUC20 Antibody, Biotin conjugated

Product Code	CSB-PA822693LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N307
Immunogen	Recombinant Human Mucin-20 protein (569-657AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Mucin-20 (MUC-20), MUC20, KIAA1359
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MUC20
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.