



NDUFAF2 Antibody, HRP conjugated

Product Code	CSB-PA836643LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N183
Immunogen	Recombinant Human NADH dehydrogenase [ubiquinone] 1 alpha subcomplex assembly factor 2 protein (117-169AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NADH dehydrogenase [ubiquinone] 1 alpha subcomplex assembly factor 2 (B17.2-like) (B17.2L) (Mimitin) (Myc-induced mitochondrial protein) (MMTN) (NDUFA12-like protein), NDUFAF2, NDUFA12L
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	NDUFAF2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.