



KIF20B Antibody, Biotin conjugated

Product Code	CSB-PA839399LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96Q89
Immunogen	Recombinant Human Kinesin-like protein KIF20B protein (1607-1731AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinesin-like protein KIF20B (Cancer/testis antigen 90) (CT90) (Kinesin family member 20B) (Kinesin-related motor interacting with PIN1) (M-phase phosphoprotein 1) (MPP1), KIF20B, KRMP1 MPHOSPH1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	KIF20B