



KIF20B Antibody, Biotin conjugated

Product Code	CSB-PA839399LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q96Q89
Immunogen	Recombinant Human Kinesin-like protein KIF20B protein (1607-1731AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinesin-like protein KIF20B (Cancer/testis antigen 90) (CT90) (Kinesin family member 20B) (Kinesin-related motor interacting with PIN1) (M-phase phosphoprotein 1) (MPP1), KIF20B, KRMP1 MPHOSPH1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	KIF20B
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.