



CD99L2 Antibody, HRP conjugated

Product Code	CSB-PA840569LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TCZ2
Immunogen	Recombinant Human CD99 antigen-like protein 2 protein (26-145AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CD99 antigen-like protein 2 (MIC2-like protein 1) (CD antigen CD99), CD99L2, MIC2L1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD99L2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.