



CD99L2 Antibody, Biotin conjugated

Product Code	CSB-PA840569LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8TCZ2
Immunogen	Recombinant Human CD99 antigen-like protein 2 protein (26-145AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	CD99 antigen-like protein 2 (MIC2-like protein 1) (CD antigen CD99), CD99L2, MIC2L1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	CD99L2