



NLRP7 Antibody, FITC conjugated

Product Code	CSB-PA840998LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WX94
Immunogen	Recombinant Human NACHT, LRR and PYD domains-containing protein 7 protein (89-187AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NACHT, LRR and PYD domains-containing protein 7 (Nucleotide-binding oligomerization domain protein 12) (PYRIN-containing APAF1-like protein 3), NLRP7, NALP7 NOD12 PYPAF3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	NLRP7