



NLRP7 Antibody, Biotin conjugated

Product Code	CSB-PA840998LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8WX94
Immunogen	Recombinant Human NACHT, LRR and PYD domains-containing protein 7 protein (89-187AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	NACHT, LRR and PYD domains-containing protein 7 (Nucleotide-binding oligomerization domain protein 12) (PYRIN-containing APAF1-like protein 3), NLRP7, NALP7 NOD12 PYPAF3
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	NLRP7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.